

Purpose

This Crash Pad serves as a soft-landing zone for children who seek sensory input or engage in therapeutic movement activities. It also functions as a large, supportive bean bag for reading, resting, or quiet time. Commonly used in schools, clinics, and supervised home environments, the pad promotes safe exploration, sensory regulation, and relaxation.

Proper Use

- Intended for body-first landings, specifically on the trunk area (torso, hips, or side).
 - Diving, flipping, or falling from floor level is acceptable when landing on the trunk, not on feet or hands.
 - Use only under adult or professional supervision.
 - Place on a flat, stable surface, free from surrounding hazards.
 - Allow the crash pad to fully expand (up to 48 hours after opening) before use.
-

Warnings & Restrictions

- Adult supervision is required at all times, especially for individuals with limited motor control or increased injury risk.
 - **Do NOT** jump or land feet-first or hands-first, this may cause injury due to foam displacement.
 - **Do NOT** dive or flip from elevated surfaces (e.g., furniture, platforms, or ladders); foam displacement may reduce impact absorption, increasing the risk of injury.
 - The TamperShield zipper helps deter access but is not childproof, keep it securely closed at all times. Children must never climb inside.
 - Internal foam pieces are not intended for removal, play, or ingestion.
 - Ensure the crash pad is placed on a non-slip surface to prevent shifting during use. For added safety, use on carpeted floors or consider applying a non-skid material to the bottom of the pad.
 - This product is not certified for high-impact activities or professional gymnastics use.
-

Maintenance & Safety Checks

- Regularly inspect for wear, tears, or damage to the cover or seams.
- Discontinue use if foam fails to return to shape or if the cover is compromised.
- Store in a dry, indoor space, away from sharp objects and direct sunlight.