

# **APPLICATION & FITTING GUIDE**

Convoluted Foam



The Skil-Care Foam Pressure Relieving Heel Protector eliminates pressure by off-loading the heel for the prevention and treatment of heel ulcers. Two interiors are available: smooth and convoluted foam. Customized for a unique fit, each heel protector comes with an auxiliary pad. The heel protectors are manufactured of latex-free polyurethane foam and have a heel insert that is covered with a low friction fabric.



## Application of Foam Heel Protector

#### STEP 1

Open up the heel protector (Fig. 1). Set the auxiliary pad aside for future use. It can be used later for customization.



STEP 2

Place patient's foot inside the heel protector. Make certain that the heel is suspended over the air cavity with toes facing upward **(Fig. 2).** 



Note: Foam Heel Protector is currently grey in color
Customizing the Fit



LONGER LEGS

Place the auxiliary pad slightly outward from heel protector to comfortably accommodate a longer leg. Cut the adhesive backing so that the adhesive is in contact with the heel protector.



Secure the auxiliary pad in a horizontal position outside the heel protector. Then secure the straps (Fig. 3).



FOOT DROP

Place the auxiliary pad in a vertical position behind the sole of the foot. Note, ensure that the heel is suspended over the air cavity.

# Care Instructions for the Foam Heel Protector

· Hand Washing: use disinfectant or soap and water.

• Machine Wash and Dry: place in a mesh laundry bag and do not dry above 180°F (82°C). After washing and drying an antiseptic spray can be used.

Item No.	Description	Size	U/M
503400	Foam Heel Protector - Smooth	Universal	Each
503450	Foam Heel Protector - Convoluted	Universal	Each
503098	Foot Stabilizer - Foam	Universal	Each

## STEP 3

Thread the Velcro<sup>®</sup> straps through the D-rings and secure the closures. (Fig. 3). Velcro<sup>®</sup> closures should be snug but not tight.



#### STEP 4

To ensure a proper fit confirm that the heel is over the heel suspension opening and the straps are comfortably tightened (Fig. 4).



## Caution \_

- 1. Use on bedridden patients only.
- 2. Do not allow patient to walk or stand on the floor while wearing the protector.
- 3. Have a medical professional remove heel protector periodically to inspect the skin.
- 4. Do not over-tighten the straps.
- 5. Assure that the foot is placed in the protector so that the heel is off-loaded.
- When used with a Deep Vein Thrombosis (DVT) device, assure that the tubing is not in contact with the skin and the air flow is not compromised.

### \*Meets CA 117 flammability standards



503098 Foot Stabilizer (Sold Separately)

