

Re-Order #	Description	Size	U/M
751570	Inflatable SquarePod Air Cushion w/LSII	16"W x 16"D x 2.5"H	Each
751571	Inflatable SquarePod Air Cushion w/LSII	18"W x 16"D x 2.5"H	Each
751572	Inflatable SquarePod Air Cushion w/LSII	18"W x 18"D x 2.5"H	Each
751573	Inflatable SquarePod Air Cushion w/LSII	20"W x 16"D x 2.5"H	Each
751580	Replacement Pump	N/A	Each



Purpose



This adjustable air cushion is designed with 25 inflatable square air pods that evenly distribute weight, enhancing circulation and reducing pressure points for greater comfort. It arrives partially inflated and can be easily customized to your desired firmness using the included air pump. The air cushion is designed to support a weight capacity of up to 500 lbs. The cushion’s low-shear, water-resistant flexible cover helps it conform to your body, providing effective relief from fatigue, numbness, and discomfort during prolonged periods of sitting—ideal for use in clinical settings, rehabilitation environments, home or while traveling.



Easy To Use

Inflation

1. Open cover zipper to the point labeled “Air Valve” and expose the air valve.
2. Open the air valve cap. (Figure A)
3. Insert air valve into air pump tube. (Figure B)
4. Inflate cushion by drawing air pump cylinder back (Figure C) and then pushing it forward (Figure D). The amount of air used to inflate the cushion will vary with the pressure relief and positioning needs of the individual user. It is up to the healthcare professional to determine the proper cushion inflation for each user.
5. Remove the air valve from the air pump tube and close valve cap (Figure E).



Deflation

1. Follow steps 1 and 2 above.
2. Squeezing the air valve to deflate cushion as long as valve is open.
3. Close valve by stop squeezing, close the cap (Figure E).

Caution

When cleaning do not use concentrated chlorine, instead use soapy water or diluted chlorine. The Cushion Air should never be machine washed or dried or placed in a sterilizer or autoclave.

Storage

WARNING: These Cushions should not be stored where the temperature is above 35°C, 95°F for any prolonged time (over 30 days is prolonged). They should not be left in a very hot vehicle for any length time exceeding 2 hours or in direct sunlight on very hot days for any length of time. It is recommended, but not required that these cushions be stored in a temperature/humidity controlled environment, and in a partially inflated state.

Cover and Cushion Care Instructions

Cover

The cover is water resistant and could be wiped clean with mild disinfectants. If necessary, the cover can be laundered at temperatures below 180°F (82°C).

Cushion

Once the cover is removed, the cushion should be hand washed with soap and water, hand washing or laundry detergent. To fully disinfect the product, wipe it clean with a diluted bleach solution (approximately one part to 9 parts of water) or follow the manufactures directions on diluting the solution. After the disinfectant is applied, wait for approximately 10 minutes and then wash with soap and water again.