

REF	Size	Waist Measurements	Unit
911450	X-Small	26" 28" 66" - 71cm	Each
911452	Small	28" 30" 71" - 76cm	Each
911454	Medium	30" 34" 76" - 86cm	Each
911456	Large	34" 38" 86" - 96cm	Each
911458	X-Large	38" 42" 96" - 106cm	Each
911460	XX-Large	42" 46" 106" - 116cm	Each
911470	Replacement Vinyl Covered Foam Pads		Pair



### Purpose

Skil-Care's Hip-Ease features high-density foam pads that are positioned over the hip to reduce impact in the event of a fall.

### Easy to Use

1. Unfasten all hook and loop straps.
2. Place waistband around users waist. The label should be positioned at the back and on the inside of the waistband.
3. Secure hook and loop closures at the front of the waist. The fit should be snug enough to prevent the device from shifting, but not so tight as to be uncomfortable.
4. Position the side panels containing the foam pads over the left and right hips. The pads are placed directly over the hip joints.
5. Fasten the hook and loop leg straps. These should be secure enough to prevent the panels from shifting, but not so tight as to be uncomfortable.



### Caution

- Skil-Care Corporation makes no guarantee, express or implied, that wearing Hip-Ease will protect the wearer against hip trauma.
- Skin under the waist and leg bands should be assessed on a regular basis.
- Hip-Ease should be changed and laundered after each incontinence episode.
- The foam pads are sealed in vinyl pouches to protect them against water or urine absorption. If the pouch seal is broken or the vinyl torn, the pads will absorb water or urine. The water or urine soaked foam may promote bacteria growth. Test the integrity of each pouch. If you hear or feel liquid or air escaping, the pouch is damaged. Discard damaged pouch and replace with reorder number 911470. Remove pad from hip-ease before laundering clean pads with mild disinfectant, or PDI wipes.

### Easy to Care

- Skil-Care Hip-Ease can be laundered and dried at temperatures up to 180°F (82°C). Remove pads before laundering, clean pads w/mild disinfectant, or PDI wipes.
- Secure hook and loop closures prior to laundering to prevent lint from adhering to hook material. If closures do not secure due to lint, remove lint from hook material using stiff brush.

