




Texture Panel (A)



Texture Panel (B)



Texture Panel (C)

Re-Order Number	Description	Size	Units
912480	Sensory Squares Sheet with Texture Panel (A)	40"W x 41"H Sheet, 12" x 12"	Set
912484	Sensory Squares - Sheet Only	40"W x 41"H	Each
912481	Sensory Squares Textures, Set of 9 pcs - Panel (A)	12" x 12"	Each
912482	Sensory Squares Textures, Set of 9 pcs - Panel (B)		Each
912483	Sensory Squares Textures, Set of 9 pcs - Panel (C)		Each

Purpose

The sensory squares provide a useful tool to help the person with autism or neurodevelopmental disorders to evaluate various textures, colors and materials. The Sensory Squares consist of nine separate 12" x 12" squares (panel A) each one varies in color and texture and a 40" x 41" backdrop sheet. Two additional sets of different squares are available; each set consists of nine separate squares (panel B & C) that are different in texture and color than panel (A). Using the system, the caregiver can design several useful tasks such as testing their skill to generalize and associate. One example is, how many blue, yellow and white squares are present?




Fig.1

Easy to Use

The display sheet can be used on a table, floor or affixed to the wall with hook, nails or screws. Orient the sheet with the hook and loop strips horizontally. To apply the squares to the sheet, have the user carefully match the hook and loop on the squares with that on the sheet (Fig.1). After the squares are placed, the user can be tested on their association and discrimination skills (i.e. Number of smooth, rough, color squares etc.). If the user has a hypersensitivity to any of the textures the caregiver could introduce the surfaces in limited doses. Additional square set "B" and "C" offer different textures, colors and materials and will offer many other task opportunities.

Easy Care

The nylon sheet can be laundered the vinyl squares can be wipe-cleaned with an anti-bacterial solution and the remaining squares can also be laundered. The nylon sheet and non-vinyl squares can be dried at temperatures that do not exceed 180°F (84°C). Do not use bleach.

Caution

Although all of the materials are child friendly, they should not be placed in the mouth. We also suggest that this product be used under the supervision of an adult.