


Bed Rail Wedge Pads



Reorder#	Description	Length	Unit
<input type="checkbox"/> 401235	Full-Rail Wedge Pads	70"	Pair
<input type="checkbox"/> 401240	Half-Rail Wedge Pads	35"	Pair


Purpose

Wedge design fills in the gap between the mattress and side rail to prevent entrapment of arms and legs. Made with extra-thick foam, these pads provide added cushioning for seizure-risk patients.

Easy to Use

1. Place pad on patient side of bed rail. Rest bottom of wedge on mattress. Position the extended section of the pad in the gap between the mattress and bed rail. See illustration above.

2. Bring straps around the outside of the rail, secure buckles, and pull loose ends of straps to take up slack. See illustration at right.



Easy Care

Skil-Care's Bed Rail Wedge Pads wipe clean with a damp cloth. A mild spray-type cleanser may be used. Do not use a cleanser that contains bleach.