

Laundered Products

Removable Fabrics Include:

- (LP1) * **Low Shear I** - A two-way stretch, breathable tricot fabric that protects against skin damage caused by friction and shear force. Fluid-proof, it has a laminated polyurethane backing.
- (LP2) * **Polyester** - A knitted fabric that stretches to reduce skin-damaging shear forces.
- (LP3) * **Cozy Cloth** - An ultra soft polyester microfiber that wicks away moisture, reduces friction and shear forces.
- (LP4) * **Nylon** - Incontinent-proof, durable and laminated with a polyurethane backing.
- (LP5) * **Polyester Mesh** - A knit, open weave polyester that is cool to the skin.
- (LP6) * **Super Soft Loop Material** -
- (LP7) * **Synthetic Sheepskin** - A fluff polyester that simulates the characteristics of sheepskin. The product helps reduce pressure sores and allows air flow between the fiber.

Cleaning:

- * The fabrics above can be laundered at a temperatures below 180°F (82°C).
- * Use lower temperature during the dry cycle to extend product life.



Wipe Cleanable Products

Fabrics Include:

- (WC1) * **Low Shear II** - A two-way stretch, breathable polyurethane fabric with a laminated polyester coating that reduces skin-damaging shear forces.
- (WC2) * **Vinyl (UnSupported)** - A heavy duty vinyl that is used for all Skil-Care welded (heat-sealed) products. The material is waterproof and anti-microbial.
- (WC3) * **Supported Vinyl (Sewn)** - A vinyl coated fabric that offers superior strength, durability, and cleanability for sewn products.
- (WC4) * **Neoprene** - A soft, pliable rubber composite with a knit fabric backing that provides maximum skin protection from abrasions and bruising.
- (WC5) * **Non-Skid Fabric** - A fabric coated with a non-skid urethane.
- (WC6) * **EVA Foam (EZ Landing Fall Mat)** - A one inch molded floor mat comprised of shock absorbing EVA foam.
- (WC7) * **Fiber Glass (Smokers Apron)** - A silicon coated woven fiber glass that is fire proof.
- (WC8) * **Non-Skid Matting** - A polyester woven fabric with a dipped rubber coating.
- (WC9) * **Nylon** - Incontinent-proof, durable and laminated with a polyurethane backing.

Cleaning:

- * Wipe fabrics with either soap and water, liquid disinfectant or disinfectant wipes.
- * Rinse thoroughly to remove any residual chemicals or soap.
- * Make sure fabric is completely dry before using.
- * Disinfect blood contamination with a 1:10 dilution of household bleach (5.25% sodium hypochlorite) as recommended by the CDC (Center of Disease Control, US Department of Health and Human Services, February 1989).

Cautions:

- * DO NOT use cleaners that contain Phenol or Benzyl.
- * DO NOT use iodophor type disinfectants, such as Betadine, because they will stain the fabrics listed above.
- * DO NOT use household bleach or concentrated cleaning agents, disinfectants or germicides without diluting the product according to the manufacturer's instructions.
- * DO NOT allow any cleaning agent to remain in contact with the fabrics for an extended time

Foam Only Products

Fabrics Include:

- (FO1) * **Vel Foam** - Disposable only
- (FO2) * **UnCovered Foam** - At the users risk, foam products sometimes can be laundered at low temperatures (160°F, 71°C) when placed in a mesh laundry bag.

Cleaning:

- * At the users risk, foam products sometimes can be laundered at low temperatures (160°F, 71°C) when placed in a mesh laundry bag.

Disposal: If the fabric has been contaminated, properly dispose of per facility policy for BIOHAZARDOUS material