

Introduction

The nursing home dilemma: Provide optimum safety conditions for the resident while adhering to government regulations which may be well intended but sometimes not beneficial. Specifically, in most situations, bed-rails are considered to be restraints.

Alarms have been proven to be effective tools for the nursing staff but 90db buzzers are sometimes disturbing to the residents. For this reason, the Skil-Care MultiPro and BedPro alarm systems alert the caregiver either through a nurses call connection or an audio alert with a pleasant music sound. In addition to the bi-level buzzer. The ChairPro alarm was designed only with a buzzer for those who don't want the music option.

Skil-CareTM
CORPORATION

29 Wells Ave.
Yonkers, NY 10701



ChairPro™ ALARMS

- Auto-Reset
- Indicator Light
- Battery Test
- Low Battery Indicator
- Battery Included
- High/Low Volume
- Buzzer Tone Only
- Cable Strain Relief
- Bumpers
- Mountable
- One Year Guarantee

The ChairPro™ alarm systems are sold as a set which includes: Sensor, Electronic Unit, Batteries and Cable Accessories.

The MultiPro and BedPro alarm units have the same features as the ChairPro alarm unit...PLUS



MultiPro™ ALARMS

- Auto Reset
- Nurse Call Capabilities
- Indicator Light
- Battery Test
- Low Battery Indicator
- Battery Included
- High-Low Volume
- Two Buzzers/
- Two Music Tones
- Cable Strain Relief
- Bumpers
- Mountable
- One Year Guarantee

*The MultiPro™ alarm is sold as a unit and does not include any of the Skil-Care sensors.

The BedPro alarm units have the same features as the ChairPro alarm unit...PLUS



BedPro™ ALARMS

- Auto Reset
- Nurse Call Cable
- In-Room On-Off
- Indicator Light
- Battery Test
- Low Battery Indicator
- AC Adaptor
- Battery Included
- Nurse Call Capabilities
- High-Low Volume
- Two Buzzer Tones/
- Two Music Tones
- Cable Strain Relief
- One Year Guarantee
- Bumpers
- Mountable
- Y-Connector

The BedPro alarm units have the same features as the ChairPro alarm unit...PLUS

Magnetic Alarms & Alarm Accessories




(A) Econo Alarm



(B) Personal Alarm



(C) Door/Window Alarm



(D) Stop Strip Alarm



(E) Replacement Magnetic/Cord



(F) Alarm Bag



(G) Universal Wall Bracket



(H) Nurses Call Accessory Kit



(I) Alarm Mounting Bracket and Belt



Skil-CareTM
CORPORATION

Pro Series Alarms

EASY AS 1-2-3



1

Plug It In



2

Place It Down



3

Test It Out

Skil-CareTM
CORPORATION

29 Wells Ave, Yonkers, NY 10701
office: 914.963.2040, toll-free: 800.431.2972

Simple...

Skil-Care alarms were designed to be used by nurses, assistants and aids with minimum training. Auto-Reset replaces ON/OFF switches, battery compartment are easily accessible, and all functions are clearly identified



ChairPro Alarm

USE: For any chair or bed

Options: Auto Reset, Indicator Light, Battery Test, Low Battery Indicator, High-Low Buzzer Tone, Cable Strain Relief, Bumpers, Screw or Strap/Velcro®, Mountable Bracket, One Year Guarantee, Free Battery Included.

Note that this unit can be mounted to any surface (Chair, Bed or Wall) with Skil-Care's optional brackets and straps.

Reliable...

Skil-Care's alarms are manufactured using precision components, state of the art SMT technology and reliability testing. Our alarms are incased in a durable ABS housing, protected by shock absorbing bumpers resulting in a dependable warning system.



MultiPro Alarm

USE: For any chair or bed

Options: Auto Reset, Indicator Light, Battery Test, Low Battery Indicator, High-Low Buzzer Tone, Cable Strain Relief, Bumpers, Screw or Strap/Velcro®, Mountable Bracket, One Year Guarantee, Free Battery Included.

PLUS Nurse Call Capabilities and Two Music Tones

Note that this unit can be mounted to any surface (Chair, Bed or Wall) with Skil-Care's optional brackets and straps.

Support...

Although Skil-Care alarms are simple to use, some training may be necessary. Among alarm manufacturers, Skil-Care has one of the largest field representative organizations and technical support is available throughout the United States.



BedPro Alarm

USE: For any chair or bed

Options: Auto Reset, Indicator Light, Battery Test, Low Battery Indicator, High-Low Buzzer Tone, Cable Strain Relief, Bumpers, Screw or Strap/Velcro®, Mountable Bracket, One Year Guarantee, Free Battery Included.

PLUS Nurse Call Cable, Y-Connector, In-Room On-Off, AC Adaptor and Two Music Tones, Nurse Call Capabilities

Note that this unit can be mounted to any surface (Chair, Bed or Wall) with Skil-Care's optional brackets and straps.



(A) UnderSeat Alarm System



(B) Sensor Pad Alarm System & Pads



(C) EZ Release Belt Alarm System



(A) UnderMattress Sensor Pads
(B) Early Warning Sensor Pads



(C) OverMattress Sensor Pad



(D) Leveling Pad Alarm



(A) UnderMattress Alarm System



(B) OverMattress Alarm System



(C) Early Warning Alarm System



(D) Seat Belt Alarm System



(E) Lift-Off Lap Cushion Alarm System



(F) Gel-Foam Cushion Alarm System



(E) UnderSeat Sensor



(F) Chair Sensor Pad




(G) EZ Release Belt Sensor
(H) Seat Belt w/Buckle Sensor




(D) Two Year Sensor Strip



(E) Auxiliary Connector Cable



(F) Bed Alarm Bracket




(G) Gel-Foam Toilet Seat Cushion Alarm Pad



(H) Leveling Pad Alarm



(I) Clip-On Alarm Bracket



(I) Nurse Call Accessory Connector Kit



(J) Lift-Off Lap Cushion Sensor



(K) Gel-Foam Cushion and Sensor

Skil-Care CORPORATION 29 Wells Ave, Yonkers, NY 10701
office: 914.963.2040, toll-free: 800.431.2972

Our sales force is available for individual product demonstrations and inservicing. For any questions regarding a distributor or sales representative please contact our customer service department or visit us online @ www.skil-care.com

If you would like to be removed from our mailing list please contact us at customerservice@skil-care.com