


SKIL-FOAM BED WEDGES

Re-Order #	Cover/No Cover	Size			Unit
		H	D	W	
<input type="checkbox"/> 553010	Cloth cover	7.5"	27"	24"	Pair
<input type="checkbox"/> 553020	No cover	7.5"	27"	24"	Pair
<input type="checkbox"/> 553030	Cloth cover	10"	27"	24"	Pair
<input type="checkbox"/> 553040	No cover	10"	27"	24"	Pair
<input type="checkbox"/> 553050	Cloth cover	12"	27"	24"	Pair
<input type="checkbox"/> 553060	No cover	12"	27"	24"	Pair



These multi-purpose foam wedges comfortably elevate your patient's head or feet. More effective than pillows or rolled towels, they resist slipping and sliding for secure positioning. Sold with or without covers, Skil-Foam Wedges are available in three sizes.

COVER LAUNDERING INSTRUCTIONS

The cover may be laundered and dried at temperatures up to 180°F (82°C). Bleaching is not recommended, Bleach will, over time, destroy the thread and fibers and shorten the useful life of the cover.

Skil-Care Corporation . . . When Safety Matters and Comfort Counts.