

Contact our customer service department for our 2010 Product Catalog.
See our full line of nursing and therapy products.





SKIL-Care
CORPORATION
from
New Products
2010

Advanced design for nursing and therapy

Visit us online at www.skil-care.com • 29 Wells Ave, Yonkers, NY 10701 • 800.431.2972



EZ Landing
Fall Mat



Cushions

Gel • Foam • Air

Stability Plus Gel-Foam Cushion

This is the ultimate cushion that combines stability, comfort pressure relief and positioning

- Separate chambers to assure stability and a raisedommel to control sliding
- Water-based gel to provide pressure relief, equalization and therapeutic cooling
- High-density foam base to prevent bottoming out
- Low-Shear I launderable outer cover reduces skin damaging friction


ReOrder #	Description	Size	U/M
751046	Stability Plus Cushion w/LSI Cover	16" x 16" x 2.5"	Ea.
751047	Stability Plus Cushion w/LSI Cover	18" x 16" x 2.5"	Ea.
751048	Stability Plus Cushion w/LSI Cover	20" x 16" x 2.5"	Ea.



Anti-Thrust Visco Foam Cushion

- Visco elastic foam molds to anatomy providing maximum pressure relief
- Leg troughs promote leg alignment for improved posture
- Raised front controls sliding

ReOrder #	Description	Size	U/M
753160	Anti-Thrust Visco Foam Cushion	18" x 16" x 2"	Ea.



Pressure Check Light/Heavy Foam Cushion

- **Dual-sided** cushion provides foam cells that adjust to the resident's anatomy
- Side with larger pods designed for heavier residents (150lbs or more)
- Side with smaller pods designed for lighter residents (149lbs or less)
- Both sides provide maximum weight distribution for pressure relief
- Low-Shear II wipe-clean outer cover reduces skin damaging friction

ReOrder #	Description	Size	Guarantee	U/M
753184	Pressure-Check Light/Heavy Cushion	18" x 16" x 3"	6 months	Ea.
753186	Pressure-Check Light/Heavy Cushion w/LSII Cover	18" x 16" x 3"	6 months	Ea.



Air Lift Seat Cushion

Pre-inflated economic cushion

- Provides decubitus prevention and comfort
- Pressure relieving air floatation surface
- Pre-inflated, easy to inflate or deflate
- Wipe-clean vinyl
- Meets flammability requirements of CA 117



Bariatric Available

ReOrder #	Description	ReOrder #	Description	ReOrder #	Description	Size	Color	U/M
756210	Cushion Only	756215	Cushion w/LSII Cover	756219	Cushion Cover w/Straps	17" x 17"	Blue	Ea.
756225	Cushion Only	756230	Cushion w/LSII Cover	756234	Cushion Cover w/Straps	19" x 19"	Grey	Ea.
756235	Cushion Only	756237	Cushion w/LSII Cover	756239	Cushion Cover w/Straps	22" x 28"	Green	Ea.
751580	Optional hand pump available for easy inflation (purchased separately)							Ea.

Foam Air Cushion w/Coccyx Cutout

Air chamber seating surface is easily adjusted to user's needs

- Inflates or deflates to meet resident's pressure reduction/comfort requirements
- Low-Shear II wipe-clean outer cover reduces skin damaging friction

ReOrder #	Description	Size	U/M
914247	Foam Air Cushion w/Coccyx Cutout	16" x 16" x 2.5"	Ea.
914248	Foam Air Cushion w/Coccyx Cutout	18" x 16" x 2.5"	Ea.
914249	Foam Air Cushion w/Coccyx Cutout	20" x 16" x 2.5"	Ea.



Hand Pump Included

Skil-Care Covers:

Low-Shear I: Two way stretch, Pique Knit, absorbent top fabric, laminated to a polyurethane water proof under layer, anti-bacterial, launderable cover.

Low-Shear II: Two way stretch, top coated polyurethane, over a polyester knit, water proof anti-bacterial, wipe-clean surface.

Wheelchair Accessories

Upper Body

Wheelchair Flip Tray

E-Z exiting arm tray

- Wheelchair tray flips out of the way for easy exiting
- Hard, soil and scratch-resistant plastic
- Side and armrest padded when in the down or flipped position
- Two cup holders for drinks, eye glasses or small items

ReOrder #	Description	Size	U/M
705031	Wheelchair Flip Tray - Right	23" X 19.5" X .5"	Ea.
705032	Wheelchair Flip Tray - Left	23" X 19.5" X .5"	Ea.



Backrest Seat Cushion

- Provides simultaneous pressure relief to the back and seating area
- Includes X-Gel cushion in seating area for pressure relief, cooling/comfort
- Low-Shear I launderable outer cover reduces skin damaging friction

ReOrder #	Description	Size	U/M
914366	Wheelchair Backrest Seat Combo w/X-Gel Cushion	16" x 16"	Ea.
914368	Wheelchair Backrest Seat Combo w/X-Gel Cushion	18" x 16"	Ea.
914443	Wheelchair Backrest Seat Combo w/X-Gel Coccyx Cushion	16" x 16"	Ea.
914444	Wheelchair Backrest Seat Combo w/X-Gel Coccyx Cushion	18" x 16"	Ea.
914372	Wheelchair Backrest w/pocket for 16" Cushion (cushion sold separately)	16" x 16"	Ea.
914374	Wheelchair Backrest w/pocket for 18" Cushion (cushion sold separately)	18" x 16"	Ea.



Fall Safety


Fall Mats



The EZ Landing fall mat stands out due to its easy to care for, one piece elastomeric foam surface, beveled low profile edges and imbedded nano silver that protects the mat from mold, mildew, stains and odor causing bacteria.

- Reduces impact energy and protects resident from serious injury
- Trip resistant beveled edges
- Nano-silver throughout which protects mat from stains, odor-causing bacteria, mold and mildew
- Textured top surface, no-slip, stable surface

ReOrder #	Description	Size	U/M
911585	EZ Landing Fall Mat	1" x 24" x 68"	Ea.
911590	EZ Landing Fall Mat	1" x 30" x 68"	Ea.



Beveled Edge

Bed Positioning & Safety

Lower Body • Wall Protection

Foam Leg/Thigh Elevator

Aligns the legs while simultaneously relaxing the spine

- Off-loads the heels for prevention of pressure sores under the knees
- Aligns the leg and thighs
- Covered with a soft launderable Cozy Cloth fabric

ReOrder #	Description	Size	U/M
555040	Foam Leg Elevator w/Scoop Cutouts - Small/Medium	17" x 7"	Ea.
555041	Foam Thigh Elevator w/Scoop Cutouts - Medium/Large	20" x 7"	Ea.



30-Degree Positioning Wedge w/Convoluted Foam Surface

Provides optimum positioning for pressure relief

- Provides optimal 30-degree angle of incline for relieving pressure on hipbone and coccyx
- Wedge is comfortable, high-density foam
- Low-Shear II wipe-clean outer cover reduces skin damaging friction

ReOrder #	Description	Size	U/M
914451	30° Bed Wedge w/Convoluted surface w/LSII Cover	7" x 12" x 24"	Pair
914452	30° Bed Wedge w/Convoluted surface w/LSII Cover	7" x 12" x 34"	Pair

Bed/Wall Protector

Protects the wall from costly repairs resulting from the bed bumping the wall

- Easy to clean nylon cover
- Non-slip backing
- Fits all beds from 36"-42"W

ReOrder #	Description	Size	U/M
401165	Bed/Wall Protector	3" x 6" x 38"	Ea.



Skin Protection

Heel

Super Soft Heel Protector

This superior heel protector combines the heel off loading characteristic in addition to providing comfort and control of foot drop

- Universally sized heel protector designed for treatment and prevention of heel pressure sore
- Open heel design effectively eliminates heel pressure by suspending the heel in the air and also promotes faster healing

ReOrder #	Description	Size	U/M
503410	Heel Protector	Universal	Ea.
503420	Heel Protector w/Anti-Rotation Wedge	Universal	Set
503430	Bilateral Anti-Rotation Wedge	Universal	Ea.



Our sales force is available for individual product demonstrations and inservicing. For any questions regarding a distributor or sales representative please contact our customer service department or visit us online @ www.skil-care.com
If you would like to be removed from our mailing list please contact us at customerservice@skil-care.com