



ReOrder #	Description	Size	Units
401166	Foam Mattress Extender	35"W x 6"H x 4"T	Ea.
401167		42"W x 6"H x 4"T	Ea.

Purpose

The Skil-Care mattress extenders can be placed between the mattress and Foot/Head board to help reduce dangerous gaps. Foam mattress extension fills the gap when using an 80" mattress on an 84" bed. Made of high-density foam. Durable, easy-to-clean vinyl cover.

Easy To Use

Place the foam mattress extender between the mattress and headboard or footboard.

Easy Care

The vinyl cover can be cleaned with a mild disinfectant.