

Re-Order #	Size	Units
501020	30" x 24"	Dozen
501050	30" x 40"	Dozen
501090	30" x 60"	Each
501210	30" x 30"	Roll/30ft



Purpose

Super-soft and over an inch thick, Comforpile bed pads aid in reducing the pressure and friction that contribute to d.u. formation. They also promote air circulation under your patient.

Easy Care

Comforpile Synthetic Sheepskin Pads may be washed and dried at temperatures up to 180°F (82°C). Bleaching is not recommended.

