





Re-Order #	Description	Size	Units
503100	Foot Elevator - Standard	4" Width	Each
503140	Foot Elevator - Standard	6" Width	Each
503101	Foot Elevator w/Synthetic Sheepskin Interior	4" Width	Each
503141	Foot Elevator w/Synthetic Sheepskin Interior	6" Width	Each

Purpose

Designed for effective decubitus ulcer prevention and management. The Skil-Care Foot Elevator raises your patient's heel off the bed to eliminate pressure and also provides the necessary air circulation for accelerated healing of existing decubiti.

Easy to Use

Wrap Foot Elevator around patient's ankle and adjust to a snug, comfortable fit. Secure hook-and-loop closure. There's no need to wrap a liner around the patient's ankle; the soft flannelette material will not irritate the patient's skin.

Easy Care

This device may be laundered and dried at temperatures up to 180° F (82° C). Bleaching is not recommended.

