



SkiL-Care[™]
CORPORATION

29 Wells Avenue
Yonkers, NY 10701

914-963-2040
800-431-2972

HEEL/ELBOW PROTECTORS

Re-Order #	Size	Unit
<input type="checkbox"/> 503203	Small	Pair
<input type="checkbox"/> 503204	Medium	Pair
<input type="checkbox"/> 503205	Large	Pair
<input type="checkbox"/> 503213	Small	36 Pr./Cs.
<input type="checkbox"/> 503214	Medium	36 Pr./Cs.
<input type="checkbox"/> 503215	Large	36 Pr./Cs.

Prevents friction burns, pressure sores and abrasions. The soft polyfoam pad is specially contoured so that it comfortably covers your patients' heels, elbows and knees. The oversized stockinette sleeve positions the pad and keeps it securely where it's needed. Skil-Care's Heel/Elbow Protector is also effective as protective covering for bandages and dressings.

LAUNDERING The Heel/Elbow Protectors may be washed and dried at temperatures up to 180° F (82°C). Bleaching is not recommended.

SkiL-Care[™]
CORPORATION

Advanced design for nursing and therapy.[™]

Visit us online
at Skil-Care.com