



Re-Order #	Description	Unit
556038	Bed Support Bolster System with Pad	Set



### Main Uses

This full body foam bolster positioning system provides a stable platform for side lying residents. When coccyx off-loading is needed. It also helps to prevent residents from rolling off the bed. The system also prevents residents arms and legs from being entrapped between the mattress and bed.

### Major Features

- Four bolsters, each 30" long, 4" high, for off-loading positioning and prevention of bed roll off
- Covers full length (75") and width(34") of a bed
- Strategic 30° bolsters allow side lying for optimum off-loading and positioning
- Comfortable, wipe clean, launderable nylon covering
- Easily attaches to the bed with eight quick release straps
- Foam wedges are removable for laundering of nylon bed system

### Caution

If the bed rails are lowered and this system is provided for bed roll off the resident should be inspected frequently to inspect for resident positioning.

