



Re-Order #	Description	Unit
556038	Bed Support Bolster System with Pad	Set



Main Uses

This full body foam bolster positioning system provides a stable platform for side lying residents. When coccyx off-loading is needed. It also helps to prevent residents from rolling off the bed. The system also prevents residents arms and legs from being entrapped between the mattress and bed.

Major Features

- Four bolsters, each 30" long, 4" high, for off-loading positioning and prevention of bed roll off
- Covers full length (75") and width(34") of a bed
- Strategic 30° bolsters allow side lying for optimum off-loading and positioning
- Comfortable, wipe clean, launderable nylon covering
- Easily attaches to the bed with six quick release straps
- Foam wedges are removable for laundering of nylon bed system

Caution

If the bed rails are lowered and this system is provided for bed roll off the resident should be inspected frequently to inspect for resident positioning.



Manufactured by: Skil-Care Corporation • Made in the U.S.A

