

| Re-Order Number | Units |
|-----------------|-------|
| 703121 | Each |



Purpose

Provides comfortable lateral head support. Designed to prevent hyperextension of the neck. Attaches to Skil-Care's Reclining Wheelchair Backrest and Air Flo Reclining Backrest as well as to all extended-back wheelchairs, geri-chairs, and recliners.

Easy to Use

1. Place Head Positioner on backrest so that it is just above the patient's shoulders.
2. Bring strap behind backrest and secure buckle. Tighten strap to prevent Head Positioner from slipping.

Easy Care

Skil-Care's Head Positioner has a polyester cover. This cover can be laundered at temperatures below 180°F (82°C)