



Re-Order #	Description	Size	Unit
703300	Seat & Back Air Cushion	17" x 34"	Each

Product Description

This combination seat and back air cushion is foldable and comes pre-inflated. The air flow within the cushion provides pressure relief and redistribution by cradling to the body and to the individual contours. It can be used on any seating surface and utilizes a low profile design which allows for easier seat height entry and exit. This combination air cushion can also be laid flat to provide 17"W x 34"L pressure relief surface where needed.

Application

Place the air cushion on any seating surface and secure cushion in place with tie straps included.

Easy Care

The cushions vinyl surface can be wiped clean using a damp cloth, spray-type cleanser or mild disinfectant. Do not use bleach.