



706365



703066



Foam Insert for 706365

Re-Order #	Description	Size	Unit
706365	Modular Lumbar Support	8"W x 17"L x 2.3"H	Each
703066	Modular Lumbar Support - Full Round	16"L x 5"D	Each

Purpose

The Lumbar Support is effective in maintaining proper lumbar lordosis while sitting in a chair or car. The uniqueness of this lumbar support is the pocket that houses a separate piece of foam. The lumbar support varies in thickness from 2" without the foam to 3" with the added foam insert. This feature extends the use to a wider audience. There is an adjustable quick release buckle that allows the product to be positioned in the chair.

Important Safety Precaution

It is suggested that prior to the use of any lumbar support, a medical professional be consulted for the proper size and/or application.

Easy Care

The Lumbar Support is covered in a wipe-clean fabric and a non-bleach disinfectant can be used. Do not launder.