



Re-Order	Description	Size	Units
708170	Transfer Board	30"L x 8"W x 1"H	Each

Purpose

The Skil-Care Padded Transfer Board provides safe, stable support for transfers between bed, wheelchair, sofa, or commode. It reduces fall risk and caregiver strain, making daily mobility safer and easier.

Easy to Use

Place the board securely under the user, with one end on the starting surface and the other on the destination surface. Guide the individual to slide smoothly across the board in a controlled motion. For added safety, a caregiver or therapist should assist.

Easy to Use

Wipe the surface with a standard disinfectant after each use. The smooth nylon finish makes cleaning quick and effective.