


Re-Order #	Description	Size	Unit
751010	Gel-Foam Cushion	16" x 16" x 2.5"	Each
751011	Gel-Foam Cushion	18" x 16" x 2.5"	Each
751012	Gel-Foam Cushion	20" x 16" x 2.5"	Each
751039	Thin-Line Gel Foam Cushion	16" x 16" x 1.5"	Each
751040	Thin-Line Gel Foam Cushion	18" x 16" x 1.5"	Each
751042	Thin-Line Gel Foam Cushion	20" x 16" x 1.5"	Each
751025	Gel-Foam Cushion w/Coccyx Cut-Out	18" x 16" x 2.5"	Each
754528	Gel-Foam Cushion w/Coccyx Cut-Out	18" x 16" x 3.5"	Each
751018	Gel-Foam Cushon Vinyl Inner Cover Only	16" x 16" x 2.5"	Each
751019	Gel-Foam Cushion Vinyl Inner Cover Only	18" x 16" x 2.5"	Each
751020	Gel-Foam Cushion Vinyl inner Cover Only	20" x 16" x 1.5"	Each
914481	Ultra Thin Gel-Foam Cushion w/LSI Cover	16" x 16" x 1"	Each
914482	Ultra Thin Gel-Foam Cushion w/LSI Cover	18" x 16" x 1"	Each
914483	Ultra Thin Gel-Foam Cushion w/LSI Cover	20" x 16" x 1"	Each





Purpose

Skil-Care's Gel-Foam Cushion is a combination gel and foam cushion designed to relieve pressure. The water based gel surface provides cooling action that reduces skin-damaging perspiration.

Important Safety Reminder

The Low-Shear cover has a safety strap that prevents the cushion from slipping off the seat. This strap must be secured to the wheelchair frame.

Easy Care

The vinyl cover may be wiped clean using a damp cloth and spray-type cleanser. A mild disinfectant may be used. The Low-Shear outer cover can be laundered and dried at temperatures that do not exceed 180°F (82°C). Do not use bleach.

