

**Purpose**


Skil-Care's Econo-Gel Pads are designed to relieve pressure and increase comfort for wheelchair users. High-density foam layers on top and bottom provide stability and cushioning. The gel-filled center inserts add pressure-relieving flotation.

**Important Notice**

Econo-Gel Pads have safety ties that must be secured to the wheelchair to prevent the cushion from slipping off the wheelchair seat.


**Easy Care**

The vinyl cover may be wiped clean using a damp cloth and spray-type cleanser or mild disinfectant. The polyester cover can be laundered and dried at temperatures that do not exceed 180°F (82°C). Do not use bleach.




**Econo-Gel Pads - Standard Series**

Re-Order Number	Cover	Size
751150	Vinyl w/Polyester	16" x 16" x 2"
751151	Vinyl w/Polyester	18" x 16" x 2"
751148	Vinyl w/Polyester	18" x 18" x 2"
751149	Vinyl w/Polyester	20" x 18" x 2"
751153	Vinyl w/Polyester	20" x 16" x 2"
751154	Vinyl Only	16" x 16" x 2"
751155	Vinyl Only	18" x 16" x 2"
751157	Vinyl Only	20" x 16" x 2"
751156	Vinyl Only	18" x 18" x 2"
751162	Vinyl w/Sheepskin	16" x 16" x 2"
751163	Vinyl w/Sheepskin	18" x 16" x 2"
751164	Vinyl w/Sheepskin	20" x 16" x 2"



**Econo-Gel Pads - Bimini Series**

Re-Order Number	Cover	Size
751180	Bimini Blue	16" x 16" x 2"
751181	Bimini Blue	18" x 16" x 2"



**Econo-Gel Pads - Thin-Line Series**

Re-Order Number	Cover	Size
751170	Vinyl w/Polyester	16" x 16" x 1.5"
751171	Vinyl w/Polyester	18" x 16" x 1.5"
751174	Vinyl w/Polyester	20" x 16" x 1.5"
751178	Vinyl w/Polyester	18" x 18" x 1.5"
751179	Vinyl w/Polyester	20" x 18" x 1.5"

Re-Order Number	Cover	Size
751172	Vinyl Only	16" x 16" x 1.5"
751173	Vinyl Only	18" x 16" x 1.5"
751175	Vinyl Only	20" x 16" x 1.5"
751176	Vinyl Only	18" x 18" x 1.5"
751177	Vinyl Only	20" x 18" x 1.5"