



- Superior comfort and pressure relief
- Large quantity of gel provides maximum pressure redistribution
- The gel layer is separated into two components to increase stability

- Reduces shear forces and enhances comfort
- Can be used on either side
- Wipe clean, water resistant, two-tone cover
- Cover includes tie straps that affixes cushion to wheelchair

Re-Order Number	Cover	Size
751166	LSII	16" x 16" x 2"
751167	LSII	18" x 16" x 2"
751168	LSII	20" x 16" x 2"

Purpose

The Super Soft Econo Gel Cushion is made with a high-density outer body to provide stability and comfort while the inner gel compartment provides superior pressure relief. The cushion comes with a water resistant, two-tone Low Shear stretch cover that is easy to clean. The two-tone allows the cushion to be used on either side and serves as a reminder to reposition the user.

Easy to Use

Super Soft Econo-Gel Cushion have safety ties that must be secured to the wheelchair to prevent the cushion from slipping off the wheelchair seat.

Easy Care

The stretch cover may be wiped clean using a damp cloth and spray-type cleanser or mild disinfectant.