

Skil-Care 29 Wells Avenue Yonkers, NY 10701 914-963-2040 800-431-2972

STARRY NIGHT GEL FOAM PAD

Reorder #	Size	
□751194	14" x 14" x 2½"	
□751196	16" x 16" x 2½"	
□751197	18" x 16" x 2½"	
□751198	18" x 18" x 2½"	
F	ReOrder# 751196	ReOrder# 751194

Purpose

The Starry Night Gel-Foam Pad features a surface layer of waterbased gel that relieves interface pressure while cooling the skin to reduce perspiration. The high-density foam base adds both comfort and stability.

Important notice

The Sittin' Pretty Starry-Gel Pad has safety ties that must be secured to the wheelchair to prevent the cushion (and possibly the user) from slipping off the wheelchair seat.

Easy care

Skil-Care's Sittin' Pretty Econo-Gel Pad is heat sealed to prevent urine and other fluids from seeping into the cushion. The vinyl cover may be cleaned with a spray-type cleanser. A mild dininfectant may be used.

Manufactured by: Skil-Care Corporation • Made in the U.S.A



