

Contour Cushion

Low-shear cover prevents skin-damaging friction. Incontinent-proof. Removes for easy laundering.

Safety straps secure cushion to wheelchair.

Elevated sides improve lateral stability.

Visco-elastic foam pad provides support and pressure relief beneath the coccyx.

Ample high density foam relieves pressure and increases comfort.

Raised front controls forward thrust and maintains pelvic neutrality to reduce sacral sitting problems.

Available gel pod provides additional pressure relief in the critical coccyx area.

Reorder #	Description	Size	Unit
751610	Contour Cushion w/Visco Foam Pad and LSII Cover	16" x 16" x 4"	Each
751615	Contour Cushion w/Visco Foam Pad and LSII Cover	18" x 16" x 4"	Eac
751620	Contour Cushion w/Visco Foam Pad and LSII Cover	20" x 16" x 4"	Each
751625	Contour Cushion w/Gel Pad and LSII Cover	16" x 16" x 4"	Each
751630	Contour Cushion w/Gel Pad and LSII Cover	18" x 16" x 4"	Each
751635	Contour Cushion w/Gel Pad and LSII Cover	20" x 16" x 4"	Each
751636	Contour Cushion w/Gel Pad and LSII Cover	22" x 18" x 4"	Each
751637	Contour Cushion w/Gel Pad and LSII Cover	24" x 18" x 4"	Each

Purpose

Designed for positioning and pressure relief. The raised front controls forward thrust to maintain pelvic neutrality. Elevated sides promote lateral stability. Ample foam cushioning relieves pressure and adds comfort.

Easy Care

The cloth cover may be wiped clean using a mild spray cleanser and damp cloth, or it may be removed for laundering. Wash and dry at temperatures below 180°F (82°C). Do not use bleach.