



Purpose

Our Basic Gel-Foam Cushion features a multi-layer design for all-day comfort and support. The one-inch foam top molds to your body for customized comfort, while the supportive foam base ensures stability. Inside, a dual-chamber gel pod provides exceptional pressure relief and weight distribution to reduce peak pressure points and prevent discomfort. The cushion is covered in a water-resistant, 2-way stretch Low Shear II fabric, offering skin protection by minimizing friction and easing cleaning for long-lasting use.

Benefits include:

- **Pressure Distribution:** Gel core evenly distributes weight, minimizing pressure points and promoting comfort.
- **Dynamic Support:** Dual foam layers balance soft comfort with firm support for optimal stability.
- **Enhanced Stability:** Prevents shifting to help maintain proper posture and balance.
- **Skin Protection:** Low-shear cover reduces friction and prevents skin irritation.
- **Lightweight and Easy Handling:** Vinyl-free design makes it easier to lift and reposition.

Easy to Use

Place the non-skid fabric bottom on the seating surface and secure with the straps and side-release buckle for a secure fit.

Easy Care

The cover can be wiped clean with a damp cloth and a mild cleanser or disinfecting wipes.

| | 2" Foam Base | 3" Foam Base | 4" Foam Base |
|-------------|----------------|----------------|----------------|
| Size | ReOrder Number | ReOrder Number | ReOrder Number |
| 16"W x 16"D | 752030 | 752031 | 752032 |
| 18"W x 16"D | 752033 | 752034 | 752035 |
| 20"W x 16"D | 752036 | 752037 | 752038 |
| 18"W x 18"D | 752039 | 752040 | 752041 |
| 20"W x 18"D | 752042 | 752043 | 752044 |
| 22"W x 18"D | 752045 | 752046 | 752047 |
| 24"W x 18"D | 752048 | 752049 | 752050 |