




Re-Order#	Description	Size	Unit
753150	Conform Cushion w/Low Shear II Cover	18"W x 16"D x 2"H*	Each
753155	Conform Cushion w/Coccyx Cutout w/Low Shear II Cover	18"W x 16"D x 2"H*	Each
753160	Conform Visco Anti-Thrust Cushion w/Low Shear II Cover	18"W x 16"D x 2"H*	Each
753157	Conform Wedge Cushion	18"W x 16"D x 4"D x 2"H*	Each

<sup>\*</sup>Front Low Point



## **Purpose**

Made with temperature-sensitive visco-foam, this cushion provides pressure relief, supports appropriate posture and controls forward sliding.

## Important Safety Reminder

The Conform Cushion with a cloth outer cover has straps with buckles that must be secured to the wheelchair to prevent the cushion from sliding off the seat.

## Easy Care

The fluid proof cover may be wiped clean using a damp cloth and spray-type cleanser. A mild disinfectant may be used. Do not use bleach.

This product is guaranteed for 6 months.

