

Conform Cushion
model# 753150




New Visco
Anti-Thrust Cushion
model# 753160

Re-Order #	Description	Size	Unit
753150	Conform Cushion w/Low Shear II Cover	18" x 16" x 1.5"	Each
753155	Conform Cushion w/Coccyx Cutout w/Low Shear II Cover	18" x 16" x 1.5"	Each
753160	Conform Visco Anti-Thrust Cushion w/Low Shear II Cover	18" x 16" x 2"	Each
753157	Conform Wedge Cushion	18" x 16" x 4 x 1.5"	Each

Low-Shear fluid proof cover reduces skin friction and shear and wipes clean for easy care.

Straps with buckle secure cushion to wheelchair.

Conform Cushion
w/Coccyx Cut-Out
model# 753155



Raised front controls forward sliding.

Temperature-sensitive visco-foam conforms to user for maximum weight distribution and pressure reduction.

Leg troughs promote thigh alignment for postural symmetry



Visco-Elastic Foam

Purpose

Made with temperature-sensitive visco-foam, this cushion provides pressure relief, supports appropriate posture and controls forward sliding.

Important Safety Reminder

The Conform Cushion with a cloth outer cover has straps with buckles that must be secured to the wheelchair to prevent the cushion from sliding off the seat.

Easy Care

The fluid proof cover may be wiped clean using a damp cloth and spray-type cleanser. A mild disinfectant may be used. Do not use bleach.

This product is guaranteed for 6 months.