



Re-Order #	Descriptions	Size	Units
753184	Pressure Check Light/Heavy - Foam Only	18" x 16" x 3"	Each
753186	Pressure Check Light/Heavy - LSII Cover	18" x 16" x 3"	Each

### **Purpose**

Pressure relief is essential for preventing decubitus ulcers and it is often the choice of a medical professional to choose a foam cushion to serve that purpose. With a standard flat cushion there is a build up of heat and a susceptibility to augmented shear force. A foam checkered surface addresses both issues and with the dual sized foam sections, residents of different weights could be accommodated. It is recommended that residents less than 140 pounds utilize the surface with smaller squares and alternatively the larger squares to be used for heavier residents.

### **Important notice**

Model number 753186 incorporates a durable water proof surface to protect the cushion against fluids and bacteria. If a tear or rip is noticed in the covering the cushion shouldn't be treated as fluid resistant.

### **Easy Care**

Model number 753186 can be wiped clean with standard disinfectants. Do not use rough cleaning pads on the membrane.