



- Superior pressure relief and comfort
- Resilient high density foam conforms to all body shapes
- Foam sections adjust individually to eliminate high and low spots and minimize shear force.
- Construction provides cooling channels throughout cushion
- Available with incontinent proof cover
- Cooling Channels
- Capacity: 300 lbs. 4" cushion, 250 lbs. 3" cushion

Re-Order Number	Description	Size
753178	Pressure Check Foam Seat Cushion with Waterproof Clear Film Covering	18" x 16" x 3"
753179	Pressure Check Foam Seat Cushion with Waterproof Clear Film Covering	18" x 16" x 4"
753180	Pressure Check Foam Seat Cushion with LSII Cover	18" x 16" x 3"
753181	Pressure Check Foam Seat Cushion with LSII Cover	18" x 16" x 4"
753182	Pressure Check Foam Seat Cushion with Sheepskin Cover	18" x 16" x 3"
753183	Pressure Check Foam Seat Cushion with Sheepskin Cover	18" x 16" x 4"

### Purpose

This cushion utilizes a unique cubic formation to provide pressure relief and redistribution while simultaneously minimizing shear forces.

### Easy to Use

The cushions should be placed on a chair with the cubes facing upward. For models without the cover, modifications such as a coccyx cut-out can easily be made with a knife.

### Easy Care

Clean the cushion cover with a germicidal cleaning solution or soap and water. Do not use bleach or launder. Store at room temperature.