

## Thin-Line™ ULTRA-CUSHION


Reorder #	Size
754340	18" x 16" x 2.5"
754339	16" x 16" x 2.5"

### PURPOSE

The Ultra-Cushion is recommended for wheelchair-bound individuals who require a firm-foundation cushion.

• **GEL DECREASES PRESSURE AND INCREASES COMFORT** For maximum pressure relief, Ultra-Cushions are available with a self-contained chamber of heat-dissipating gel that evenly distributes resident's weight.

• **PROVIDES COMFORT** An ample layer of pliant high-density foam offers maximum cushioning for resident comfort.



• **URINE CAN'T PENETRATE COVER** The Ultra-Cushion is protected by a factory-sealed, incontinent-proof vinyl cover.

• **ELIMINATES WHEELCHAIR SEAT HAMMOCKING** Noncompressing firm foam foundation eliminates the sling-seat hammocking that causes spinal pressure. The contoured foundation rests below wheelchair seat frame—your resident sits lower in the chair.

### IMPORTANT SAFETY FEATURE

The cloth outer cover has a safety strap that prevents the cushion from slipping off the seat. This strap must be secured to the wheelchair frame.

### COVER CARE INSTRUCTIONS

Vinyl covers may be cleaned with a mild spray-type cleaner and damp cloth or sponge. Cloth covers can be removed for laundering. Do not launder or dry cloth covers at temperatures in excess of 180°F (82°C). Bleaching is not recommended.