



**LOW-SHEAR I**



**LOW-SHEAR II**



**SYNTHETIC SHEEPSKIN**



**TWO-TONE LOW-SHEAR II**



**TWO-TONE POLYESTER**

Size	LEVEL			
	Comes with a Non-Slip Bottom			
	Low-Shear I (LSI)	Low-Shear II (LSII)	Synthetic Sheepskin	UM
16"W x 16"D x 1"- 2"H	781014	781028	781034	Ea.
16"W x 16"D x 3"- 4"H	781015	781029	781035	Ea.
18"W x 16"D x 1"- 2"H	781016	781030	781036	Ea.
18"W x 16"D x 3"- 4"H	781017	781031	781037	Ea.
20"W x 16"D x 1"- 2"H	781018	781032	781038	Ea.
20"W x 16"D x 3"- 4"H	781019	781033	781039	Ea.

Size	LEVEL		
	Two-Tone Reminder Cover		
	LSII	Polyester	UM
16"W x 16"D x 1"- 2"H	781054	781044	Ea.
16"W x 16"D x 3"- 4"H	781055	781045	Ea.
18"W x 16"D x 1"- 2"H	781056	781046	Ea.
18"W x 16"D x 3"- 4"H	781057	781047	Ea.
20"W x 16"D x 1"- 2"H	781058	781048	Ea.
20"W x 16"D x 3"- 4"H	781059	781049	Ea.

Size	WEDGED	
	Comes with a Non-Slip Bottom	
	Low-Shear II (LSII)	UM
16"W x 16"D x 4"- 2"H	781040	Ea.
16"W x 16"D x 6"- 3"H	781041	Ea.
18"W x 16"D x 4"- 2"H	781042	Ea.
18"W x 16"D x 6"- 3"H	781043	Ea.

**Launderable Products** - Launder at a temperatures below 180°F (82°C).

**Low Shear I** - A two-way stretch, breathable tricot fabric that protects against skin damage caused by friction and shear force. Fluid-proof, it has a laminated polyurethane backing.

**Polyester** - A knitted fabric that stretches to reduce skin-damaging shear forces.

**Synthetic Sheepskin** - A fluff polyester that simulates the characteristics of sheepskin. The product helps reduce pressure sores and allows air flow between the fiber.

**Wipe-Clean Products** - Wipe fabrics with either soap and water, liquid disinfectant or approved disinfecting wipes.

**Low Shear II** - A two-way stretch, breathable polyurethane fabric with a laminated polyester coating that reduces skin-damaging shear forces.