



| Re-Order # | Description | Size | Units |
|------------|-----------------------------------|------------------------------|-------|
| 905030 | Elevating Arm Support - Foam Only | 20"L x 6"W x 5.5" Max Height | Each |

Purpose

This soft arm support provides a trough for resting the arm in a comfortable and stable position. Fits either the left or right arm and can be placed on any leveled surface such as a table, bed or one of Skil-Care's wheelchair arm trays.

- · Recessed cavity
- Can be used to reduce upper extremity edema
- Can be used as a phlebotomy wedge during blood draws.
- Product dimensions: 20"L x 6"W x 5.5" Max Height

Easy to Use

Place the armrest positioner on any leveled surface or wheelchair arm tray and assure the arm is centered.