



Replacement Foam & Elastic Banding



914334K - Set

Re-Order	Description	Size	Units
914332K	Reusable Face Shield	Universal	12 per case
914333K	Replacement foam & elastic banding	Universal	12/PK
914334K	(2) Reusable Face Shields and (8) replacement foam and banding	Universal	Set

Purpose

Face shields are personal protective equipment devices used to protect the facial area while acting as a shield against transferring droplets from the nose or mouth. These face shields are perfect for a school setting. It doesn't block the mouth so you can see a teacher talking as well as it is more comfortable to wear during a long periods of time.

Composition

The shield is manufactured from a high quality Polyethylene Terephthalate Glycol material commonly referred to as PETG. It is non-porous and can be cleaned with anti-bacterial wipes. The elastic banding and foam support are removable and replaceable as recommended.

Warning: Do not remove the foam from the shield unless you are ready to replace the foam.

Easy Use

This face shield is not sterile, if required, it can be gas sterilized. Place the shield over the face and adjust it to maximize comfort. If desired, the elastic banding and foam can be peeled off after each use.

Cleaning Instructions

This face shield can be cleaned with anti-bacterial wipes. The elastic banding and foam support are removable and replaceable as recommended.

Caution

This product is not sterilized.