





Foam Support

Elastic Banding
Clear Face Shield





Replacement Foam & Elastic Banding

914334 - Sei

Re-Order #	Description	Size	Units
914332	Reusable Face Shield	Universal	12 per case
914333	Replacement foam & elastic banding	Universal	24/PK
914334	(2) Reusable Face Shield w/Extra Foam (8) & Banding (8)	Universal	Set

Purpose

Face shields are personal protective equipment devices used to protect the facial area while acting as a shield against splashes in the work environment.

Composition

The shield is manufactured from a high quality Polyethylene Terephthalate Glycol material commonly referred to as PETG. It is non-porous and can be cleaned with anti-bacterial wipes. The elastic banding and foam support are removable and replaceable as recommended. **Warning: Do not remove the foam from the shield unless you are ready to replace the**

Easy Use

This face shield is not sterile, if required, it can be gas sterilized. Place the shield over the face and adjust it to maximize comfort. If desired, the elastic banding and foam can be peeled off after each use.

Cleaning Instructions

This face shield an be cleaned with anti-bacterial wipes. The elastic banding and foam support are removable and replaceable as recommended.

Caution

This product is not sterilized.



foam.