



Re-Order Number	Description	Size	Units
914471	Wheelchair Wedge Base	18" x 16" x 2"	Each
914472	Wheelchair Wedge Base	18" x 16" x 3"	Each
914473	Wheelchair Wedge Base	18" x 16" x 4"	Each

### **Purpose**

This product when placed in a chair serves a dual purpose, 1. Namely elevating the front, 2. Converting a hammocked fabric surface to a rigid base. The rigid surface helps maintain a neutral seating position while the raised front helps prevent the resident from sliding off the wheelchair. In the event the resident is seated in a reclining back chair, the elevation helps maintain a proper 90 degree hip alignment.

### **Easy to Use**

The proper elevation should be determined by a medical professional. Once the proper model is identified the wedge base should be placed on the chair should be tied around the chair

### **Easy Care**

The plastic base can be wiped clean with soap and water or a standard disinfectant. Always inspect that the tie straps are in proper working condition.

***This product is guaranteed for one year.***