

Vibrator Mitt

ReOrder # 914560



SkiL-Care's Vibe-Mitt provides a soft vibrating massage for relaxation and sensory stimulation. It offers a choice of message textures: high-pile synthetic sheepskin on one side and velvet-like Cozy Cloth on the other. The vibrator requires two AA batteries (1.5 volts each). It is pressure operated –there's no on/off switch.

Caution : Do not mix old and new batteries, Do not mix alkaline, standard (carbon zinc) or re-chargeable (nickel-cadmium) batteries

Easy Care

The Vibe-Mitt is completely launderable. Wash and dry at temperatures up to 180⁰ F (82⁰C). **Remove vibrator before laundering.**