



Re-Order #	Description	Size	Unit
915060	Budget Bariatric Foam Cushion w/LSII Cover	20"W x 16"D x 4"H	Each
915061	Budget Bariatric Foam Cushion w/LSII Cover	20"W x 18"D x 4"H	Each
915062	Budget Bariatric Foam Cushion w/LSII Cover	22"W x 16"D x 4"H	Each
915063	Budget Bariatric Foam Cushion w/LSII Cover	22"W x 18"D x 4"H	Each
915064	Budget Bariatric Foam Cushion w/LSII Cover	24"W x 18"D x 4"H	Each
915065	Budget Bariatric Foam Cushion w/LSII Cover	24"W x 20"D x 4"H	Each
915066	Budget Bariatric Foam Cushion w/LSII Cover	26"W x 18"D x 4"H	Each
915067	Budget Bariatric Foam Cushion w/LSII Cover	26"W x 20"D x 4"H	Each
915068	Budget Bariatric Foam Cushion w/LSII Cover	28"W x 18"D x 4"H	Each
915069	Budget Bariatric Foam Cushion w/LSII Cover	28"W x 20"D x 4"H	Each

Purpose

Made with pressure reducing and protective foam for providing relief to most bariatric residents. Skil-Care has two other gel-filled bariatric cushions and a three layer (Tri-Foam) bariatric cushion for extreme (400+ lbs.) bariatric residents. In the event the resident is less than 400 lbs. this economy bariatric cushion would be effective for pressure reduction and would also be significantly less expensive.

Features

The cushion is comprised of very high quality four inch foam for added comfort, pressure redistribution and reduction of pressure sores. In addition, it is covered with a two way stretch, waterproof fabric that will aid in the reduction of shear forces and protect against incontinence. In addition, there is a safety strap for attachment to the wheelchair.

Easy Care

The fluid proof cover may be wiped clean using a damp cloth and spray-type cleanser. A mild disinfectant may be used. Do not use bleach.

Caution

The budget bariatric foam cushion has straps with buckles that must be secured to the wheelchair to prevent the cushion from sliding off the seat. In addition, other Skil-Care bariatric cushions are recommended for residents under 400 lbs.