




**Reminder**

Two-color surfaces acts as a reminder to reposition resident



Re-Order	Description	Size	Units
751140	Bio-shield Econo-Gel Cushion - Yellow	16" x 16" x 2"	Each
751142	Bio-shield Econo-Gel Cushion - Green	18" x 16" x 2"	Each



**Purpose**

The Skil-Care Bio-shield cushion is designed to provide pressure relief and add comfort with the additional benefit of bacteria-fighting properties. The cushion's outer cover is made from antifungal, antibacterial wipe clean vinyl to decrease the spread of infection. The Cushion comes with a gel-filled insert surrounded by a top and bottom layer of foam. The Two-Color surfaces aids as a reminder to reposition resident.

**Important Notice**

Comes with a safety tie to secure the cushion to a wheelchair.

**Easy Care**

Clean cushion with spray and wipe cleaning solution, avoid bleach. Do not launder cushion.