



# Contour Plus<sup>+</sup> Cushion

## Contour Plus<sup>+</sup> Wedge

Reorder #	Description	Size
<input type="checkbox"/> 751675	Contour Plus+ Cushion	18" x 16" x 3"
<input type="checkbox"/> 751680	Contour Plus+ Wedge	18" x 16" x 5" x 2"

### Purpose

Contour Plus Cushions and Wedges are designed for wheelchair users who require superior wheelchair positioning and pressure relief. Both styles feature convex bases that eliminate sling-seat hammocking to reduce pressure on the spine and coccyx.



Available wedged base provides additional slide control.

### Easy Care

The Low-Shear cover may be wiped clean, or it may be removed for laundering. The cover should be washed and dried at temperatures that do not exceed 180°F ( 182°C). Do not use bleach.