



ReOrder	Description	Size	U/M
999999	Cylinder Bolster	16.5"L x 5"D	Each

Purpose

This Cylindrical Bolster is made with a high-density foam and has an anti-bacterial fluid proof, wipe clean cover. The cover material is Low Shear II which is a two-way stretch material that decreases skin shear forces. It can be placed under the knee, ankle or upper extremity to maintain the scapular plain. It can also be used between the trunk and upper extremity during band exercises.

Easy to Use

The bolster is manufactured from high-density foam and the cover is wipe clean. Most alternate bolsters are made from a porous foam or alternatively cloth covers and rolled up towels harbor bacteria.

Easy Care

The antibacterial fabric can be easily wiped cleaned with disinfectant wipes/ liquids or soap and water. Do not launder!