

ULTRA-WEDGE



Skil-Care™

29 Wells Avenue
Yonkers, NY 10701

914 963-2040
800 431-2972

Description	Cover	Size (inches)	Reorder #	Size (inches)	Reorder #	Size (inches)	Reorder #
Foam Pad	Vinyl Only	16 x 18 x 3½ x 1½	754110 ☐	16 x 16 x 3½ x 1½	754112 ☐	16 x 20 x 3½ x 1½	754114 ☐
Foam Pad	Vinyl w/Cloth	16 x 18 x 3½ x 1½	754120 ☐	16 x 16 x 3½ x 1½	754122 ☐	16 x 20 x 3½ x 1½	754124 ☐
Gel-Foam	Vinyl w/Cloth	16 x 18 x 4 x 1½	754130 ☐	16 x 16 x 4 x 1½	754132 ☐	16 x 20 x 4 x 1½	754134 ☐
Foam Pad	Vinyl only	16 x 18 x 4½ x 1½	754160 ☐	16 x 16 x 4½ x 1½	754162 ☐	16 x 20 x 4½ x 1½	754164 ☐
Foam Pad	Vinyl w/Cloth	16 x 18 x 4½ x 1½	754170 ☐	16 x 16 x 4½ x 1½	754172 ☐	16 x 20 x 4½ x 1½	754174 ☐
Gel-Foam	Vinyl w/Cloth	16 x 18 x 5 x 1½	754180 ☐	16 x 16 x 5 x 1½	754182 ☐	16 x 20 x 5 x 1½	754184 ☐

PURPOSE

The Ultra-Wedge is recommended for wheelchair-bound individuals who require the slide prevention and positioning capabilities of a wedge cushion and the thigh rotation control and spinal pressure relief of a firm seating foundation.

- TIE STRAPS FOR RESIDENT SAFETY** Ultra-Wedges have ties to prevent them from slipping off wheelchairs.
- PRESSURE-REDUCING GEL** For maximum pressure relief, choose the Gel-Foam Ultra-Wedge. A separate chamber of heat-dissipating gel evenly distributes resident's weight.
- URINE CAN'T PENETRATE COVER** The Ultra-Wedge is protected by a heat-sealed, incontinent-proof vinyl cover. Other wedges have stitched, zippered, or envelope covers that allow urine to penetrate to the foam.
- PROVIDES RESIDENT COMFORT** An ample layer of pliant high-density foam offers maximum cushioning for resident comfort.
- ELIMINATES WHEELCHAIR SEAT HAMMOCKING** A noncompressing firm foam foundation eliminates the sling-seat hammocking that causes spinal pressure and resident sliding. Separate wheelchair leveling pads are unnecessary.
- FIXED WEDGE SHAPE** The firm foundation resists compression and maintains a constant slide-controlling angle of elevation.

IMPORTANT SAFETY FEATURE

All Ultra-Wedge covers—vinyl or cloth—have tie cords to prevent the wedge from slipping off the wheelchair seat. For resident safety, these cords must be secured to the wheelchair frame.

COVER CARE INSTRUCTIONS

Vinyl covers may be cleaned with a mild spray-type cleaner and damp cloth or sponge. Cloth covers can be removed for laundering. Do not launder or dry cloth covers at temperatures in excess of 180°F (82°C). Bleaching is not recommended.