

# CONVOLUTED FOAM PAD

Re-Order #	Cover	Size	Re-Order #	Cover	Size
☐910110	No Cover	18" x 16" x 4"	☐910112	LSII	18" x 16" x 4"
☐910115	No Cover	18" x 16" x 3"	☐910117	LSII	18" x 16" x 3"

## Purpose

Skil-Care's Convoluted Foam Pad is recommended for use on wheelchairs and geri-chairs as well as standard chairs. The pad reduces pressure, increases comfort, and improves air circulation.

## WARNINGS

- Do not expose pad to open flames, burning cigarettes, and other sources of ignition. If ignited, pad will burn rapidly and give off toxic fumes.
- This pad is intended for single patient use. Do not wash. Laundering will cause the pad to lose its fire retardant characteristics.

Skil-Care's Convoluted Foam Pad meets California 117B flammability standards.