


# COMFORT FOAM WHEELCHAIR CUSHION

Reorder #	Cover	Cushion Size (inches)
<input type="checkbox"/> 707052	Vinyl	18" x 16" x 6"
<input type="checkbox"/> 753210	LSII	16" x 18" x 2"
<input type="checkbox"/> 753230	Vinyl	16" x 18" x 2"
<input type="checkbox"/> 753232	Vinyl	16" x 16" x 2"
<input type="checkbox"/> 753240	No Cover	16" x 18" x 2"
<input type="checkbox"/> 753250	LSII	16" x 18" x 3"
<input type="checkbox"/> 753270	Vinyl	16" x 18" x 3"
<input type="checkbox"/> 753280	No Cover	16" x 18" x 3"
<input type="checkbox"/> 753285	LSII	16" x 18" x 4"
<input type="checkbox"/> 753290	Vinyl	16" x 18" x 4"
<input type="checkbox"/> 753295	No Cover	16" x 18" x 4"
<input type="checkbox"/> 753268	Vinyl	16" x 16" x 2"
<input type="checkbox"/> 753272	Vinyl	16" x 16" x 3"
<input type="checkbox"/> 753271	Vinyl	18" x 18" x 3"
<input type="checkbox"/> 753251	LSII	18" x 18" x 3"
<input type="checkbox"/>		

Maximum cushioning with minimum compression make Skil-Care High Density Foam Cushions your best choice for wheelchair and geri-chair comfort.



## COVER CARE

Nylon and vinyl covers may be cleaned with a spray-and-wipe cleaner and damp cloth. Laundering the nylon cover will cause it to lose its water repellent quality. If the nylon cover is laundered, do not wash or dry it at temperatures over 180°F (82°C). Bleaching is not recommended.