Skil-Care
 29 Wells Avenue
 914-963-2040

 Yonkers, NY 10701
 800-431-2972

SOFT BOOT

Re-Orde	er#	Size	Unit
□50316	0	Universal	Pair

Skil-Care's Soft Boots offer extended decubitus protection that provides pressure relief for your patient's heel, ankle, and lower calf. Made from plush synthetic sheepskin, these boots are comfortably secured with hook-and-loop closures. Large vent holes promote air circulation to reduce perspiration and moisture build-up.

LAUNDERING INSTRUCTIONS

The Soft Boots may be laudered and dried at temperatures up to 180° F (82° C). Bleaching is not recommended. To prolong the useful life of the hook-and-loop closures, secure them before laundering. This will prevent the hook material from picking up lint during the washing and drying process.

Advanced design for nursing and therapy.™





