


Vinyl Bed Rail Pads

Thru-View Vinyl Bed Rail Pad

Full-View Vinyl Bed Rail Pads


Full Rail (pair)



NEW Full-View Vinyl Bed Rail Pads

Reorder #	Product Description	Size
<input type="checkbox"/> 401010	Vinyl Bed Rail Pads	15" x 1" x 60"
<input type="checkbox"/> 401015	Thru-View Vinyl Bed Rail Pads	15" x 1" x 60"
<input type="checkbox"/> 401040	Vinyl Bed Rail Pads	15" x 1" x 72"
<input type="checkbox"/> 401045	Thru-View Bed Rail Pads	15" x 1" x 72"
<input type="checkbox"/> 401050	Vinyl Bed Rail Pads	15" x 1" x 80"
<input type="checkbox"/> 401100	Vinyl Bed Rail Pads	15" x 2" x 60"
<input type="checkbox"/> 401140	Vinyl Bed Rail Pads	15" x 2" x 72"
<input type="checkbox"/> 401150	Vinyl Bed Rail Pads	15" x 2" x 80"
<input type="checkbox"/> 911525	Vinyl Bed Rail Pads (Special Order)	15" x 1" x 60"
<input type="checkbox"/> 912200	Full-View Vinyl Bed Rail Pads	15" x 1" x 60"
<input type="checkbox"/> 912205	Full-View Vinyl Bed Rail Pads	15" x 1" x 72"
<input type="checkbox"/> 912210	Full-View Vinyl Bed Rail Pads	15" x 1" x 80"
<input type="checkbox"/> 401330	Thin-Line Bed Rails Pads	66" x 11.5" x .5"

Half Rail (4 Units/Package)



NEW Full-View Bed Rail Pads.

Reorder #	Product Description	Size
<input type="checkbox"/> 401080	Vinyl Bed Rail Pads (Half Rail)	9" x 1" x 28"
<input type="checkbox"/> 912220	Full-View Vinyl Bed Rail Pads (Half Rail)	15" x 1" x 29"
<input type="checkbox"/> 912225	Full-View Vinyl Bed Rail Pads (Half Rail)	15" x 1" x 40"

Purpose

Skil-Care's bed rail pads are designed to protect patients against side rail entrapment and against skin tears and bruises resulting from contact with metal rails. These pads are not intended as substitutes for proper supervision of patients or residents for whom bed rails are being used.

Application

Place bed rail pads on patient side of rail and secure Velcro[®] over rails as shown in illustrations.

Easy Care

All Skil-Care vinyl bed rail pads may be cleaned using spray-and-wipe cleanser and a damp cloth. Mild disinfectants may be used.

Note

Bed rail pad covers are non-allergenic, mildew resistant, bacteriostatic, and conform to flammability standards set forth in California Bulletin 117B.