



WHEELCHAIR ARMREST CUSHION

Reorder #	Description	Unit
□703015	Full-Arm (15")	Pair
□703016	Half-Arm (11")	Pair
□703130	Extended Full Arm (22")	Pair
□703131	Extended Half Arm (18")	Pair



PURPOSE

Foam-padded and nylon-covered, Skil-Care Armrest Cushions are designed to provide comfort and protect sensitive arms and elbows against contact with hard armrests and metal wheelchair frames.

APPLICATION

Wrap Armrest Cushions around wheelchair armrests and secure Velcro underneath armrest.

LAUNDERING

Armrest Cushions can be laundered and dried at temperatures up to 180°F (82°C). Bleaching is not recommended.